

BWTL-BALIWALA AND HOMI TRIAL LENS SET

WITH WOODEN BOX-RED



Baliwala and Homi-Trial Lens Set Model-BWTL
With Wooden Box-Red

- A.** Reduced Aperture Lenses Fitted in Metal Mounts with Powder Coated Red-Black Color
- B.** SPHERES
Convex and Concave
0.12
0.25 to 4.00 in 0.25 steps
4.50 to 6.00 in 0.50 steps
7.00 to 14.00 in 1.00 steps
16.00 to 20.00 in 2.00 steps
- C.** CYLINDERS
Convex and Concave
0.12
0.25 to 3.50 in 0.25steps.
4.00 to 6.00 in 0.50steps
- D.** 10 Prisms
1/2, 1, 2, 3, 4, 5, 6, 8, 10, 12
- F.** Accessories
Red Glass, Green Glass, Pin Hole, Slit & 2 Blank Discs
- G.** Lens Diameter-20mm Inner, 38mm Outer

Product Video:

<https://youtu.be/qPvJ4HcUGY>

Manufactured & Marketed By:

MEX
ENTERPRISES

We Share in Healthcare...
ISO 13485 :2016 & ISO 9001 :2015 Certified Company

HEAD OFFICE

MEX ENTERPRISES
ISO 13485 :2016 & ISO 9001 :2015 Certified Company
Plot No. 50/1B-50/2B, Rajharsh Colony, Near CJ Square Shopping Mall, Kolar Road, Bhopal, M.P
India-462042

FACTORY ADDRESS

MEX ENTERPRISES
ISO 13485 :2016 & ISO 9001 :2015 Certified Company
Plot No.-23B, Badiyakhedi, Industrial Park Phase-1, Badiyakhedi, Sehore, MP
India-466001

DELHI SALES OFFICE

MEX ENTERPRISES
ISO 13485 :2016 & ISO 9001 :2015 Certified Company
Building Name - Classic's, T-1895/1, Upper Ridge Road, Karol Bagh, New Delhi India- 110005

RAIPUR SALES OFFICE

MEX ENTERPRISES
ISO 13485 :2016 & ISO 9001 :2015 Certified Company
House No.-10/806, ShivanandNagar, Jhanda Chowk, Raipur Chhattisgarh India-492008

Ph.: +91 755 4256391, +91 9302778762, +91 9111429222, Email: info@mediequipbazar.com, export@mediequipbazar.com
Website: www.mediequipbazar.com

www.mediequipbazar.com