


# LM-1000

## MANUAL LENS METER



## FEATURES



## TECHNICAL SPECIFICATION

Type	Ocular, External, Reading
Target	Movable cross corona target
Vertex P, Scale	0 to +/-24 diopter marked in 0.25 steps upto 10D, in 0.50 steps
Cylinder Axis	Beyond 10D
Prism	0 to 180 Degrees in 1 Degree increment
Diopter	5 Prism Diopter in 1 Diopter increment
Lens Size	Upto 90mm
Tilting	Possible continuously variable from 30 to 90°
Finish	Staple Gray
Power Supply	3 Volt DC(Two Pencil Cell)/220V AC
Spare	2 Pencil cell dust cover 2 green bulb



+91 9302 778 762  
+91 9111 429 222



export@mediequipbazar.com  
info@mediequipbazar.com



+91 9302 778 762  
+91 9111 429 222



## MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi  
Industrial Park Phase I,  
Badiyakhedi, Sehore, M.P.,  
India - 466 001



### MANUFACTURING UNIT

Plot No.-40, KI Ansal  
Green, Kolar Road  
Bhopal-462042



### MANUFACTURING UNIT

Plot No. - 79, Commercial  
Shade, Near D-Mart,  
Bhopal, 462042



[www.mediequipbazar.com](http://www.mediequipbazar.com)



### REGISTERED OFFICE

#### MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,  
Near CI Square Shopping Mall, Kolar Road,  
Bhopal, M.P.-462042

### SALES OFFICE

#### MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,  
Upper Ridge Road, Karol Bagh,  
New Delhi - 110005

### SALES OFFICE

#### MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,  
Jhanda Chowk, Raipur CG-492008