

MAT-100

APPLANATION TONOMETER



FEATURES

Reliable and repeatable aplanation tonometry, the Medix Aplanation Tonometer is a perfect example of Indian manufacturing experience delivering great quality and value. We have applied 21st Century design techniques to an old established principle to deliver a premium quality instrument.

ACCURATE MEASUREMENT

With years of tests and clinical applications, Goldmann has been proved to be the gold standard for tonometry measurement. Based on the proven evidence, Medix MAT tonometer applied Goldmann principle in the design and manufacture. With precision manufacturing process, Medix aplanation tonometer offers accurate and consistent measurement which give a clear picture of patients' IOP situation to the doctor.

EXCELLENT MECHANICS

Full attentions are paid to every details in the manufacture of Medix tonometers.

PREMIUM OPTICS

The prism on the tonometer is the key to accurate measurement. In Medix, our skill full workers make these prism with great care and follow strict quality standards. We offer 3 pcs of prism in standard packing.



- **Measuring range: 0 - 1064 Kpa(0-8 mmHg).**
- **Halo Displacement: 1.53 x 2 =3.06 mm.**
- **Prism head diameter: 7 mm.**
- **Prism head range of activities: 3 mm.**

TECHNICAL SPECIFICATION

Range of Measurement	0mmHg~80mmHg (0 kPa~1064 kPa)
Moving Range of the Rings	1.53x2=3.06mm
Prism Diameter	7mm
Moving Range of Prism	+ 1.5mm
Dimension	290mm x 170mm x 110mm
Net Weight	500g (T Type) 800g (R Type)



+91 9302 778 762
+91 9111 429 222



export@mediequipbazar.com
info@mediequipbazar.com



+91 9302 778 762
+91 9111 429 222



MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi
Industrial Park Phase I,
Badiyakhedi, Sehore, M.P.,
India - 466 001



MANUFACTURING UNIT

Plot No.-40, KI Ansal
Green, Kolar Road
Bhopal-462042



MANUFACTURING UNIT

Plot No. - 79, Commercial
Shade, Near D-Mart,
Bhopal, 462042



www.mediequipbazar.com



REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
Near CI Square Shopping Mall, Kolar Road,
Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
Upper Ridge Road, Karol Bagh,
New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,
Jhanda Chowk, Raipur CG-492008