

MEDIX
MEDICAL

MEX
ENTERPRISES
We Share in Healthcare...

MD- 03 & 05

THREE & FIVE STEP SLIT LAMPS

WITH HALOGEN / LED ILLUMINATION



www.mediequipbazar.com

MD-03 HS & MD-05 HS

The MEDIX slit lamp is a stereoscopic microscope that emits a focused beam of light with variable height, width, and angle.

MEDIX MD-03 HS (H)/(L), MD-05 HS (H)/(L), delivers ultimate value with a well-designed 3/5 step optical system, provide clearest view, crisp and details of anterior and posterior segment of patient's eye.

All parts in contact with patient are Antitoxic and Anti-irritation and an efficient halogen/led illumination.

TECHNICAL SPECIFICATION

SPECIFICATIONS-BINOCULAR MICROSCOPE

Viewing Oculars	: Binocular, Galilean Type
Magnification Change	: Five/Three Steps Drum Rotation Type Magnification
Eye Pieces	: Wide Field 12.5X
Objective	: 0.4x, 0.6x, 1x, 1.6x, 2.5x/ 0.6x, 1x, 1.6x
Total magnification	: 5x, 7.5x, 12.5x, 20x, 30x/ 7.5x, 12.5x, 20x

ILLUMINATION UNIT

Slit Image Rotation	: 0 to 180 degree
Tilting Illumination	: 5,10, 20 degree
Filter Disc	: Cobalt Filter, Red Free, Natural Density, milky white and Open Aperture Diaphragm
Slit Disphragm Disc	: Six Apertures of 12, 8, 5, 3, 1, 0, 0.2mm
Halogen Lamp	: 12 Volt 50 Watts/5 Watts Led

STANDARD ACCESSORIES

Replacing Mirror	: Plastic Breath Shield
Testing Rod	: Spare Lamp (Halogen Lamp)/Led-No Extra Led



+91 9302 778 762
+91 9111 429 222



export@mediequipbazar.com
info@mediequipbazar.com



+91 9302 778 762
+91 9111 429 222



MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi
Industrial Park Phase I,
Badiyakhedi, Sehore, M.P.,
India - 466 001

REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
Near CI Square Shopping Mall, Kolar Road,
Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
Upper Ridge Road, Karol Bagh,
New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,
Jhanda Chowk, Raipur CG-492008



MANUFACTURING UNIT

Plot No.-40, KI Ansal
Green, Kolar Road
Bhopal-462042



MANUFACTURING UNIT

Plot No. - 79, Commercial
Shade, Near D-Mart,
Bhopal, 462042



www.mediequipbazar.com



www.mediequipbazar.com