

# SRS-400

## STREAK RETINOSCOPE



## FEATURES

- Enhanced light output makes refraction easier.
- Superior streak quality that enables easier and quicker refractions.
- Smooth, non-abrasive knurled finished handle provide better grip.
- Available in Rechargeable and Dry cell options.
- Dedicated carry pouch makes it ultra portable.

## TECHNICAL SPECIFICATION

### SREAK RETINOSCOPE

<b>Type:</b>	Streak Retinoscope
<b>Illumination:</b>	3.5V Halogen Illumination
<b>Power Supply:</b>	Rechargeable: 3.7V Rechargeable Battery with Charger Dry Cell: 2X Dry Cell 1.5V AA Size
<b>Handle:</b>	Knurled Metal Handle
<b>Head:</b>	Durable Plastic Head
<b>Carry Case:</b>	Dedicated carry case with packing

---



+91 9302 778 762  
+91 9111 429 222



export@mediequipbazar.com  
info@mediequipbazar.com



+91 9302 778 762  
+91 9111 429 222



## MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi  
Industrial Park Phase I,  
Badiyakhedi, Sehore, M.P.,  
India - 466 001



### MANUFACTURING UNIT

Plot No.-40, KI Ansal  
Green, Kolar Road  
Bhopal-462042



### MANUFACTURING UNIT

Plot No. - 79, Commercial  
Shade, Near D-Mart,  
Bhopal, 462042



[www.mediequipbazar.com](http://www.mediequipbazar.com)



### REGISTERED OFFICE

#### MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,  
Near CI Square Shopping Mall, Kolar Road,  
Bhopal, M.P.-462042

### SALES OFFICE

#### MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,  
Upper Ridge Road, Karol Bagh,  
New Delhi - 110005

### SALES OFFICE

#### MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,  
Jhanda Chowk, Raipur CG-492008