


MP-400 MANUAL PHOROPTER



TECHNICAL SPECIFICATION

Sphere Power	Range	+16.75 to -19.00 diopters
	Step	0.25 diopters
Cylinder Power	Range	0.00 to -6.00 diopters(0.00 to -8.00 diopters with additional lens of -2.00 diopters)
	Step	0.25 diopters
Cylinder Axis	0 to 180°, 5°step	
Cross Cylinder	±0.25 diopters	
Rotary Prisms	0 to 20, 1 step	
Auxiliary Lens Dial	(O)-Open aperture (two positions)	
	(R)-Retinoscopic lens, +1.50D	
	(P)-Polarizing lens 45-left eye, 135-right eye	
	(WMV) or (RMV)-Maddox rod, vertical white-left eye red-right	
	(WMV) or (RMV)-Maddox rod, horizontal white-left eye red-right	
	(RL)-Red lens	
	(+.12)-+0.12D sphere	
	(PH)-Pin hole	
	(10 L) or (6 U)-10 base-in-left eye	
	6 base-up-right eye	
	(Dissociating prisms)(±.50)-±0.50 fixed cross cylinder	
(OC)-Occluder		
PD Adjustment	50~75MM, 1MM steps	
Standard Accessory Lenses	wo-.12D and -2.00 cylinder lenses, two plano lenses	
Standard Accessories	1 near point card with holder and reading rod, dust cover, face shield, accessory case for lenses	
Dimensions	31.8*29.3*9.6(CM)	
Net Weight	4.5KGS	



+91 9302 778 762
+91 9111 429 222



export@mediequipbazar.com
info@mediequipbazar.com



+91 9302 778 762
+91 9111 429 222



MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi
Industrial Park Phase I,
Badiyakhedi, Sehore, M.P.,
India - 466 001



MANUFACTURING UNIT

Plot No.-40, KI Ansal
Green, Kolar Road
Bhopal-462042



MANUFACTURING UNIT

Plot No. - 79, Commercial
Shade, Near D-Mart,
Bhopal, 462042



www.mediequipbazar.com



REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
Near CI Square Shopping Mall, Kolar Road,
Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
Upper Ridge Road, Karol Bagh,
New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,
Jhanda Chowk, Raipur CG-492008