

MEDIX
MEDICAL

MEX
ENTERPRISES
We Share in Healthcare...

PORTONO+

NON CONTACT TONOMETER

HAND-HELD NCT



www.mediequipbazar.com

14
mmHg

Portono+

PORTONO+ HAND HELD NCT

Non-Contact Tonometer without anesthesia, pain and contamination concerns:

Unlike conventional contact-type tonometer, Portono+ is safe to use without the need for anesthesia, pain or contamination

Self-measuring: IOP can be measured at home

The Portono+ is light, easy to carry self-measurement at anytime, anywhere use

Easy-to-use Operation

The Portono+ is a convenient device for any hospital need, anyone can use it easily with very little practice.

Patient Comfort


The highly portable unit allows measurements to be taken wherever the patient is most comfortable; whatever environment they are in. The sleek product profile means the clinician can maintain visual contact with the patient at all times, enabling them to reassure the patient throughout the measurement process

CORE TECHNOLOGY

- Pneumatic Nozzle design
- IOP reading and Distance calculation using multiple Infra-Red Leds
- Left/Right Alignment identification using Multiple Irs
- Central alignment guide light
- User friendly OLED display screen
- Easy to use compact design

MODE SELECT

MEX
Portono+



TECHNICAL SPECIFICATION

Measurement Mode:	AUTO
Measurement Range:	5-50 mmHg
Measurement Value:	1 mmHg step
Measurement Value Storage:	20 (Twenty) Times
Power Saving Function:	Turns off in 5 Min no operation Power button return
Display:	OLED Panel 1.54 inch
Power:	12V DC 1A Input
Dimension:	177 (W) x 54 (D) x 110 (H) mm
Weight:	440 gms
Contains:	Portono+ Console, Power Supply, User Manual and Packaging

MEDIX
Portono+



+91 9302 778 762
+91 9111 429 222



export@mediequipbazar.com
info@mediequipbazar.com



+91 9302 778 762
+91 9111 429 222



MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi
Industrial Park Phase I,
Badiyakhedi, Sehore, M.P.,
India - 466 001

REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
Near CI Square Shopping Mall, Kolar Road,
Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
Upper Ridge Road, Karol Bagh,
New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,
Jhanda Chowk, Raipur CG-492008



MANUFACTURING UNIT

Plot No.-40, KI Ansal
Green, Kolar Road
Bhopal-462042



MANUFACTURING UNIT

Plot No. - 79, Commercial
Shade, Near D-Mart,
Bhopal, 462042



www.mediequipbazar.com



www.mediequipbazar.com