


OT-FM

OPHTHALMIC OPERATION TABLE



OT-FM Motorized Operation Table

- The table is designed and made specially according to ophthalmologists' needs and can apply to various operations for ophthalmology.
- Tabletop lift adopts electric motorized type, easy to operate using foot switch, silent movement of motor.
- All accessories (IV pole, tray, drape holder) made of 304 stainless steel, easy to clean.
- Nice appearance and stable performance can be regarded as an ideal ophthalmic operating table.
- Adjustable headrest height, surgeon arm rest height and bed tilt (reverse Trendelenburg system)
- Equipped with manual surgeon chair, makes the doctor more comfortable and convenient in operation.




Motorized Up-Down



Stable Base with Mobility



Bed Tilt System



Adjustable Head Rest



304 SS Accessories



Easy Foot Switch Control



TECHNICAL SPECIFICATION

OT TABLE

Maximum Height:	860 mm
Minimum Height:	660 mm
Vertical Adjustment Range:	200 mm
Total Length:	1900 mm
Total Width:	600 mm
Bed Thickness:	90 mm
Base Length:	1000 mm
Base Width:	540 mm
Trendelenburg:	10 degree
Reverse Trendelenburg:	6 degree
Power:	100~240 AC 50Hz
Auxiliary Connection:	3 Auxiliary Sockets
Weight Capacity:	200 Kilograms
Gross Weight:	98 Kilograms
Standard Accessories:	Manual Surgeon Chair, IV Pole, Instrument Tray with stand, Drape holder, fastening belts, foot switch

SURGEON CHAIR

Maximum Height:	670 mm
Minimum Height:	520 mm
Vertical Adjustment Range:	150 mm



+91 9302 778 762
+91 9111 429 222



export@mediequipbazar.com
info@mediequipbazar.com



+91 9302 778 762
+91 9111 429 222



MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi
Industrial Park Phase I,
Badiyakhedi, Sehore, M.P.,
India - 466 001



MANUFACTURING UNIT

Plot No.-40, KI Ansal
Green, Kolar Road
Bhopal-462042



MANUFACTURING UNIT

Plot No. - 79, Commercial
Shade, Near D-Mart,
Bhopal, 462042



www.mediequipbazar.com



REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
Near Cl Square Shopping Mall, Kolar Road,
Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
Upper Ridge Road, Karol Bagh,
New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,
Jhanda Chowk, Raipur CG-492008