

RU-FM

OPHTHALMIC REFRACTION UNIT



RU-FM MOTORIZED REFRACTION UNIT

- Stylish, comfortable and durable patient chair with high quality seat and upholstery.
- Motorized Up-Down chair movement function using foot switch and control panel.
- Imported DC motor for silent and maintenance free use.
- Sliding table top for two instrument provides ease of refraction process without moving the patient.
- Built-in drawer for trial set provide easy access of lenses and protects the lenses from dust.
- Built in Vision chart/drum stand, IO stand, near vision stand and torch/DO/Retino tray
- Comes with a high quality manual doctor stool of same color combo.




Patient's Chair



Manual Doctor Stool



Table Top for 2 Instrument



Control Panel



Utility Pole



Easy Foot Switch Control

OPTIONAL ACCESSORIES



MANUAL PHOROPTER



DIGITAL PHOROPTER



PHOROPTER ARM



LED VISION CHART



TRIAL LENS SET



NEAR VISION DRUM



DIRECT OPH.



INDIRECT OPH.

TECHNICAL SPECIFICATION

CHAIR UNIT

Type:	Fully Motorized
Up-Down:	Motorized
Recliner:	Motorized
Tabletop for:	Two Instruments
Chair Elevation Stroke:	200 mm
Seat Height:	560 to 760 mm
Max. Area Required:	1450x1250mm
Rotation:	Yes - 0-180 degree
Load Lifting Capacity:	135 Kgs
Auxiliary Output:	230V 50Hz * 4
Power Supply:	230V 50/60Hz 5A
Gross Weight:	240 Kgs

DOCTOR STOOL

Maximum Height:	670 mm
Minimum Height:	520 mm
Vertical Adjustment Range:	150 mm



+91 9302 778 762
+91 9111 429 222



export@mediequipbazar.com
info@mediequipbazar.com



+91 9302 778 762
+91 9111 429 222



MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi
Industrial Park Phase I,
Badiyakhedi, Sehore, M.P.,
India - 466 001



MANUFACTURING UNIT

Plot No.-40, KI Ansal
Green, Kolar Road
Bhopal-462042



MANUFACTURING UNIT

Plot No. - 79, Commercial
Shade, Near D-Mart,
Bhopal, 462042



www.mediequipbazar.com



REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
Near CI Square Shopping Mall, Kolar Road,
Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
Upper Ridge Road, Karol Bagh,
New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No. -10/806, Shivanand Nagar,
Jhanda Chowk, Raipur CG-492008