

MST-400 SCHIOTZ

APPLANATION TONOMETER



MST-400 SCHIOTZ APPLANATION TONOMETER

- Medix Schiötz tonometer model MST-400
- Professor Schiötz's Theories, best tool for measuring intraocular pressure. New Generation Indian technology made product perfect.
- Precision measurement on a scale of 0 to 20 subdivisions and 0 to -1. sub-divisions.
- Perfect reading of the scale with pointer.
- MST-400 complete set Supplied with three weight(5.5 g, 7.5 g, 10 g) & premium wood carry case



TECHNICAL SPECIFICATION

Dimensions: Height: 130 mm, Width: 45 mm

Weight: 26 g Approx

Scale: 0 to 20 graduation marks

0 to -1 graduation mark

One graduation mark corresponds a stroke of 0.05 mm

Operating Temperature: 10 degree C to + 40 degree C

Relative Humidity: 30% up to 70 % relative humidity, non-condensing

Storage Temperature: -20 degree C to +70 degree C



+91 9302 778 762
 +91 9111 429 222



export@mediequipbazar.com
 info@mediequipbazar.com



+91 9302 778 762
 +91 9111 429 222



MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi
 Industrial Park Phase I,
 Badiyakhedi, Sehore, M.P.,
 India - 466 001



MANUFACTURING UNIT

Plot No.-40, KI Ansal
 Green, Kolar Road
 Bhopal-462042



MANUFACTURING UNIT

Plot No. - 79, Commercial
 Shade, Near D-Mart,
 Bhopal, 462042



www.mediequipbazar.com



REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
 Near CI Square Shopping Mall, Kolar Road,
 Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
 Upper Ridge Road, Karol Bagh,
 New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,
 Jhanda Chowk, Raipur CG-492008