

# MEDIX TS-1100(IMPORTED) TRIAL LENS SET

WITH ALUMINUM CASE SILVER METAL RIM WITH TRIAL FRAME & CROSSCYLINDER



Medix Trial Lens Set Model-TS-1100 (Imported) with Aluminum Case Silver Metal Rim With Trial Frame and Crosscylinder.

- A. 268 pcs of Glass Lens.
- B. Lens Diameter 36.5mm(38mm outer)
- C. Spheres:40 Pairs Each of Concave and Convex-  
0.25D to 6.00D in 0.25 steps  
6.50D to 10.00D in 0.50 steps  
11.00D to 15.00D in 1.00 steps  
16.00D to 20.00D in 2.00 steps
- D. Cylinders:20 Pairs Each of Concave and Convex-  
0.25D to 4.00D in 0.25steps.  
4.50D to 6.00D in 0.50steps
- E. Prisms:12 Pieces  
0.5(2) 1.0 to 10.0 in 1.0steps
- F. Accessories:16 Pieces  
RF GF BL PL(2) MR(2) FL CL(2) PH(2)  
SS PF Crosscylinder x0.25x0.50

**Product Video:**

<https://youtu.be/qPvJ4HcUGY>

Manufactured & Marketed By:

**MEX**  
ENTERPRISES  
*We Share in Healthcare...*  
ISO 13485 :2016 & ISO 9001 :2015 Certified Company

**HEAD OFFICE**

**MEX ENTERPRISES**  
ISO 13485 :2016 & ISO 9001 :2015 Certified Company  
Plot No. 50/1B-50/2B, Rajharsh  
Colony, Near Cj Square Shopping  
Mall, Kolar Road, Bhopal, M.P  
India-462042

**FACTORY ADDRESS**

**MEX ENTERPRISES**  
ISO 13485 :2016 & ISO 9001 :2015 Certified Company  
Plot No.-23B, Badiyakhedi,  
Industrial Park Phase-1,  
Badiyakhedi, Sehore, MP  
India-466001

**DELHI SALES OFFICE**

**MEX ENTERPRISES**  
ISO 13485 :2016 & ISO 9001 :2015 Certified Company  
Building Name - Classic's, T-1895/1,  
Upper Ridge Road, Karol Bagh, New  
Delhi India- 110005

**RAIPUR SALES OFFICE**

**MEX ENTERPRISES**  
ISO 13485 :2016 & ISO 9001 :2015 Certified Company  
House No.-10/806, ShivanandNagar,  
Jhanda Chowk, Raipur Chhattisgarh  
India-492008

Ph.: +91 755 4256391, +91 9302778762, +91 9111429222, Email: [info@mediequipbazar.com](mailto:info@mediequipbazar.com), [export@mediequipbazar.com](mailto:export@mediequipbazar.com)  
Website: [www.mediequipbazar.com](http://www.mediequipbazar.com)

[www.mediequipbazar.com](http://www.mediequipbazar.com)