

# ID-05 WIRELESS

## LED INDIRECT OPHTHALMOSCOPE





The illumination rheostat is located on the headband providing easily accessible fully controllable light. It has the brightest illumination that can be turned down to as low as 5% of maximum if required.

Fully adjustable optics allow you to achieve the optimal lighting and view of your patients eyes. With the ID-05 all you need is right at your fingertips. Three spot sizes are integrated along with three filters, Yellow filter , Red free filter and cobalt blue filter. The Tiltable mirrors allows you to adjust the light beam + or - 4 degrees.

## TECHNICAL SPECIFICATION

<b>Illumination:</b>	LED
<b>Light Source:</b>	3W LED
<b>Illumination Control:</b>	ON/OFF, Linear Variation (Low to High)
<b>Filters:</b>	Heat Absorbing, Green, Cobalt Blue & Yellow
<b>Apertures:</b>	Small, Medium & Large
<b>Focus Distance:</b>	200 mm to 800 mm
<b>Field:</b>	Clear Circular
<b>Illumination Area:</b>	75mm, 55mm, 20mm (500mm Distance)
<b>Illumination Spot Size:</b>	5.5 mm, 3.5 mm, 1.5 mm
<b>Battery Backup:</b>	More than 6 Hours
<b>Battery Type:</b>	Dry Cell x 2 (AA Size)/ Rechargeable
<b>Battery Charging Time:</b>	Less Than 2 Hours
<b>Net Weight:</b>	520 Gms
<b>Carrying Case:</b>	Rigid Carrying Case







+91 9302 778 762  
 +91 9111 429 222



export@mediequipbazar.com  
 info@mediequipbazar.com



+91 9302 778 762  
 +91 9111 429 222



**MANUFACTURING UNIT**

Plot No. 23-B, Badiyakhedi  
 Industrial Park Phase I,  
 Badiyakhedi, Sehore, M.P.,  
 India - 466 001

**REGISTERED OFFICE**

**MEX ENTERPRISES**

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,  
 Near CI Square Shopping Mall, Kolar Road,  
 Bhopal, M.P.-462042

**SALES OFFICE**

**MEX ENTERPRISES**

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,  
 Upper Ridge Road, Karol Bagh,  
 New Delhi - 110005

**SALES OFFICE**

**MEX ENTERPRISES**

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,  
 Jhanda Chowk, Raipur CG-492008



**MANUFACTURING UNIT**

Plot No.-40, KI Ansal  
 Green, Kolar Road  
 Bhopal-462042



**MANUFACTURING UNIT**

Plot No. - 79, Commercial  
 Shade, Near D-Mart,  
 Bhopal, 462042



www.mediequipbazar.com

